

## **Streptavidin Fluorescent Conjugates**

| Catalog No. | Streptavidin Fluorescent Conjugate  | Concentration | Volume |
|-------------|-------------------------------------|---------------|--------|
| NEL720      | Streptavidin-Fluorescein            | 1 mg/ml       | 1 ml   |
| NEL721      | Streptavidin-Texas Red <sup>®</sup> | 1 mg/ml       | 1 ml   |
| NEL722      | Streptavidin-Coumarin               | 1 mg/ml       | 1 ml   |

### **Introduction**

Fluorescence detection is becoming one of the most powerful tools in research laboratories. As equipment (e.g. fluorescence microscopes, CCD cameras, imaging software) to measure signals becomes more sophisticated, fluorescence technology is becoming more widespread.

The streptavidin fluorescent conjugates offered by NEN Research Products have been successfully used in fluorescent immunohistochemistry and *in situ* hybridization. More specifically, these conjugates are compatible with the **TSA-Indirect** kit (NEL700), which is a signal amplification system (1-6) for immunohistochemistry. The panel of fluorescent conjugates allows for multicolor fluorescence detection. Streptavidin-fluorescein gives a greenish yellow color, streptavidin-Texas Red<sup>®</sup> gives a red color and streptavidin-coumarin gives a blue color.

The following is a table of excitation and emission wavelengths for fluorescein, Texas Red<sup>®</sup> and coumarin:

| <u>Fluorophore</u>     | <u>Excitation</u> | <u>Emission</u> |
|------------------------|-------------------|-----------------|
| fluorescein            | 494 nm            | 517 nm          |
| Texas Red <sup>®</sup> | 593 nm            | 612 nm          |
| coumarin               | 402 nm            | 443 nm          |

### **Preparation of Conjugates**

The conjugates are prepared by reacting streptavidin and the appropriate activated fluor. The fluorescent conjugate is purified by gel filtration and evaluated for fluor/streptavidin ratio. A fluor/streptavidin ratio of greater than 1 has given excellent results in immunohistochemistry applications.

### **Composition**

The streptavidin fluorescent conjugates are stored in a phosphate buffer containing BSA and sodium azide. The conjugates are at a concentration of 1 mg/ml.

### **Stability and Storage**

The products are shipped on blue ice and should be stored at 4°C upon receipt. Under the proper storage conditions the fluorescent conjugates will be stable for 6 months.

### **References**

- (1) Bobrow, MN, Harris, TD, Shaughnessy, KJ, and Litt, GJ, (1989) Catalyzed reporter deposition, a novel method of signal amplification, Application to immunoassays, *J. Immunol. Meth.* 125, 279–285.
- (2) Bobrow, MN, Shaughnessy, KJ, Litt, GJ, (1991) Catalyzed reporter deposition, a novel method of signal amplification, II. Application to membrane immunoassays, *J. Immunol. Meth.* 137, 103–112.
- (3) Adams, JC, (1992) Biotin Amplification of Biotin and Horseradish Peroxidase Signals in Histochemical Stains, *J. Histochem. Cytochem.* 40, 1457–1463.
- (4) Berghorn, KA, Bonnett, JH, Hoffman, GE, (1994) cFos Immunoreactivity is Enhanced with Biotin Amplification, *J. Histochem. Cytochem.* 42, 1635–1642.
- (5) Raap, AK, van de Corput, MPC, Verenne, RAW, van Gijlswijk, RPM, Tanke, HJ, Wiegant, J, (1995) Ultra-Sensitive FISH using Peroxidase-Mediated Deposition of Biotin-or Fluorochrome Tyramides, *Human Mol Genet* 4, 529-534.
- (6) Kerstens, HMJ, Poddighe, PJ, Hanselaar, AGJM, (1995) A Novel *in situ* Hybridization Signal Amplification Method, based on the Deposition of Biotinylated Tyramine, *J. Histochem. Cytochem.* 43/4, 347-352.

### **Hazard Information**

This product is considered to be non-hazardous. Although the product is classified as non-hazardous, we strongly recommend using prudent laboratory practices: avoiding unnecessary contact, use of gloves, eye protection, lab coats, etc. when using this or any laboratory reagent.

Texas Red<sup>®</sup> is a trademark of Molecular Probes, Inc.